

HUMAN MEDITEK

Low Temperature Hydrogen Peroxide Sterilizer



Company History Timeline

2000.03

Founded Human Meditek Co. Ltd.

2001. 12

Acquired of Manufacturing Authorization for Plasma Sterilizers from KFDA

2002.10

Designated as 'Next Generation World-Class Product' by Ministry of Commerce, Industry and Energy (MOCIE), Korea

2002.12

Grand Prize Korea Industrial Design Awards

2003.08

Acquired CE Certificate (MDD 94/42/EEC), ISO 9001, EN 46001, ISO13485

2003.10

Global Partners Signing Ceremony & Next Generation Plasma Sterilizer Official Launching Ceremony for "HMTS-80"



2005. 04

Launched HMTS-SES

2007. 06

HMTS-Series Registered in Health Canada and TGA (Australia)

2007.09

Launched "HMTS-142"















HMTS-Series registered in **53** Countries world-wide Participated in over **45** Medical Device Exhibitions and Sterilization Congresses for 15 years

AMERICAS



Canada Mexico Brazil Ecuador Chile Columbia Venezuela Peru

EUROPE



Belgium Greece Italy Poland Portugal Spain Switzerland Austria Turkey Norway Bulgaria Czech Russia Kazakhstan Sweden Netherlands Denmark Slovakia UK France

THE MIDDLE EAST



Iran Saudi Arabia UAE Syria Egypt Jordan Kuwait Lebanon

ASIA / OCEANIA



Korea Japan China Thailand India Indonesia Singapore Taiwan Vietnam Sri Lanka Philippines Australia New Zealand

2008.05

Innovative Product Award by Prime Minister of Korea

2008.05

Awarded "Export Tower Prize" from KITA

2009.02

Launched "HMTS-80E" and "HMTS-40E"

2009. 12

Acquired ISO 14973 Certificate (SMP) **2011. 02** Launched "HMTS-30E"



2013.07

Acquired Olympus Endoscope Sterilization Certificate

2013. 11

Awarded "The Majestic Five Continent Award" for Quality and Excellence

2017. 01

Developed ENDO-Mode System

2017. 10

Acquired Pentax Endoscope Sterilization Certificate

2017. 11

Acquired Prion Inactivation Certificate

2018

NEW HMTS-142 Launch Scheduled




















Low Temperature Hydrogen Peroxide Sterilizer HMTS SERIES



- Up to 89% chamber volume usage
- Up to 18 cycles per one sterilant bottle
- Safe and eco-friendly drain system
- Thousands of devices listed for compatibility
- Double door option available









HMTS-30E

HMTS-30E offers speed and versatility in a compact low temperature sterilizer. The rapid 25 minute SURFACE Mode with 30L usable chamber volume is the ideal choice for anywhere in the hospital.

Dimensions and Weight

Size Width 740 mm

Height 825 mm (With Cart) 1350 mm

Depth 725 mm

Weight 160 kg

Sterilization Chamber

Total Volume 34 L (1.2 cubic feet)
Usable Volume 30 L (1.0 cubic feet)

Dimensions (W x H x D) 340 mm x 170 mm x 600 mm

Shelf Strength 25 kg

Shelf Dimensions (W x D) 320 mm x 560 mm

STERILIZATION MODE OVERVIEW

SURFACE MODE



STANDARD MODE



ADVANCE MODE



Quick Surface Sterilization Cycle

- General metal/non-metal devices requiring surface sterilization, or sterilization of mated stainless steel or titanium surfaces
- Non lumen rigid, semi-rigid endoscopes

Conventional Flexible Lumen Sterilization Cycle with Extra Loads

- Ten(10) Polyethylene and Teflon lumen instruments with:
 Internal diameter of 1 mm or larger and a length of 500 mm or shorter
- One single channel flexible endoscope with: Internal diameter of 1mm or larger and a length of 850mm or shorter

Versatile Flexible/Rigid Sterilization Cycle with Extra Loads

• **Ten(10)** Stainless steel lumen instruments with:

Internal diameter of 1 mm or larger and a length of 500 mm or shorter

• Five(5) Stainless steel lumen instruments with:

Internal diameter of 0.77 mm or larger and a length of 500 mm or shorter

One shelf of the chamber can be used in SURFACE mode.

Max. loading weight: 1.89kg (4.1lbs)

Up to two shelves of the chamber may be used in STANDARD mode.

Max. loading weight: 3.72kg (8.2lbs)

Up to two shelves of the chamber may be used in ADVANCED mode.

Max. loading weight: 3.72kg (8.2lbs)

No.of cycle/H₂O₂ bottle 18 cycles No.of cycle/H₂O₂ bottle
12 cycles

No.of cycle/H₂O₂ bottle 9 cycles

^{*} All medical devices to be sterilized must be washed and dried completely prior to sterilization in HMTS system.



HMTS-80E

HMTS-80E is the most popular model in the HMTS lineup. With a large chamber and four cycle options to process a wide variety of medical devices, HMTS-80E enables customers to enhance overall productivity and enable a high standard of patient care.

Dimensions and Weight

Size Width 750 mm

Height* 1670 mm *Including built-in wheels

Depth 835 mm

Weight 430 kg

Sterilization Chamber

Total Volume 80 L (2.83 cubic feet)

Usable Volume 71 L (2.5 cubic feet)

Dimensions (W \times H \times D) Ø400 mm x 650 mm (D)

Shelf Strength 35 kg

Shelf Dimensions (W x D)

Upper: 380 mm x 635 mm Lower: 260 mm x 635mm

STERILIZATION MODE OVERVIEW

SURFACE MODE



ENDO MODE



STANDARD MODE



ADVANCE MODE



Quick Surface Sterilization Cycle

- General metal/non-metal, non-lumened devices requiring surface sterilization, or sterilization of mated stainless steel or titanium surfaces

Non lumen rigid, semi-rigid endoscopes

Efficient Flexible Endoscope

Sterilization Cycle

- Two(2) Single channel flexible endoscopes with: Internal diameter of 1mm or larger and a length of 850mm or shorter
- One(1) Dual channel flexible endoscope with: Internal diameter of 1mm or larger and a length for one channel of 998mm or shorter and the other channel of 850mm or shorter.

Conventional Flexible Lumen Sterilization Cycle with Extra Loads

• Ten(10) Polyethylene and Teflon lumen instruments with: Internal diameter of 1 mm or larger and a length of 500 mm or shorter

• One(1) Long Polyethylene and Teflon lumen instruments with: Internal diameter of 1mm or larger and a length of 1600mm or shorter

Versatile Flexible/Rigid Sterilization Cycle with Extra Loads

- Ten(10) Stainless steel lumen instruments with: Internal diameter of 1 mm or larger and a length of 500 mm or shorter
- Five(5) Stainless steel lumen instruments with: Internal diameter of 0.77 mm or larger and a length of 500 mm or shorter

One shelf of the chamber can be used in SURFACE mode.

Max. loading weight: 4.45kg (9.8lbs)

Up to two shelves of the chamber may be used in ENDO mode.

*Only Flexible Endoscopes with accessories may be loaded in ENDO mode Up to two shelves of the chamber may be used in STANDARD mode.

Max. loading weight: 8.9kg (19.6lbs)

Up to two shelves of the chamber may be used in ADVANCED mode.

Max. loading weight: 8.9kg (19.6lbs)

No.of cycle/H2O2 bottle 14 cycles

No.of cycle/H2O2 bottle

No.of cycle/H2O2 bottle

No.of cycle/H2O2 bottle

^{*} All medical devices to be sterilized must be washed and dried completely prior to sterilization in HMTS system.



HMTS-142

HMTS-142 is the most advanced low temperature hydrogen peroxide sterilizer from Human Meditek. With the large 142L chamber and optimized cycle options, HMTS-142 allows the customers to process more loads* and a wide variety of medical devices quickly, safely and efficiently.

 $\mbox{\ensuremath{\mbox{1}}}\xspace$ to three (3) flexible surgical endoscopes with accessories.

Dimensions and Weight

Size Width 833 mm

Height* 1688 mm *Including built-in wheels

Depth 1016 mm

Weight 550 kg

Sterilization Chamber

Total Volume

142 L (5.01 cubic feet)

Usable Volume

120 L (4.24 cubic feet)

Dimensions (W \times H \times D)

425 mm x 425 mm x 790 mm

Shelf Strength

(1 to 3 shelves)

35 kg

Shelf Dimensions (W x D)

Upper: 415 mm x 775 mm Lower: 410 mm x 775mm

STERILIZATION MODE OVERVIEW

SURFACE MODE —



ENDO MODE



STANDARD MODE



ADVANCE MODE



Quick Surface Sterilization Cycle

- General metal/non-metal, non-lumened devices requiring surface sterilization, or sterilization of mated stainless steel or titanium surfaces
- Non lumen rigid, semi-rigid endoscopes

Flexible Endoscope Sterilization Cycle

- Three(3) Single channel flexible endoscopes with: Internal diameter of 1mm or larger and a length of 850mm or shorter
- Two(2) Dual channel flexible endoscope with: Internal diameter of 1mm or larger and a length for one channel of 998mm or shorter and the other channel of 850mm or shorter

Conventional Flexible Lumen Sterilization Cycle with Extra Loads

- Ten(10) Polyethylene and Teflon lumen instruments with: Internal diameter of 1 mm or larger and a length of 1,200 mm or shorter
- One(1) Long Polyethylene and Teflon lumen instruments with: Internal diameter of 1mm or larger and a length of 4,000mm or shorter

Versatile Flexible/Rigid Sterilization Cycle with Extra Loads

- Ten(10) Stainless steel lumen instruments with: Internal diameter of 1 mm or larger and a length of 500 mm or shorter
- Five(5) Stainless steel lumen instruments with : Internal diameter of 0.77 mm or larger and a length of 500 mm or shorter

One shelf of the chamber can be used in SURFACE mode.

Max. loading weight: 7.4kg (16.3lbs)

Up to three shelves of the chamber may be used in ENDO mode.

*Only Flexible Endoscopes with accessories may be loaded in ENDO mode

Up to two shelves of the chamber may be used in STANDARD mode.

Max. loading weight: 14.8kg (32.6lbs)

Up to two shelves of the chamber may be used in ADVANCED mode.

Max. loading weight: 14.8kg (32.6lbs)

No.of cycle/H2O2 bottle
12 cycles

No.of cycle/H₂O₂ bottle

No.of cycle/H2O2 bottle
7 cycles

No.of cycle/H2O2 bottle 6 cycles

^{*} All medical devices to be sterilized must be washed and dried completely prior to sterilization in HMTS system.



HMTS-142D

HMTS-142D is a pass-through low temperature hydrogen peroxide sterilizer designed for the Central Sterile Services Departments with a need for an isolated clean room. The automatic sliding door allows for easy loading and increase productivity.

Dimensions and Weight

Size Width 1050 mm

Height* 1688 mm *Including built-in wheels

Depth 1075 mm

640 kg Weight

Sterilization Chamber

Total Volume

142 L (5.01 cubic feet)

Usable Volume

120 L (4.24 cubic feet)

Dimensions (W x H x D)

425 mm x 425 mm x 790 mm

Shelf Strength

35 kg

ADVANCE MODE

Shelf Dimensions (W x D)

Upper: 415 mm x 775 mm Lower: 410 mm x 775mm

STERILIZATION MODE OVERVIEW

SURFACE MODE

Quick Surface Sterilization Cycle

- General metal/non-metal, non-lumened devices requiring surface sterilization, or sterilization of mated stainless steel or titanium surfaces
- Non lumen rigid, semi-rigid endoscopes

ENDO MODE

Efficient Flexible Endoscope

- Three(3) Single channel flexible endoscopes with:
- endoscope with: Internal diameter of 1mm or larger and a length for one channel of 998mm or shorter and the other

Sterilization Cycle

- Internal diameter of 1mm or larger and a length of 850mm or shorter • Two(2) Dual channel flexible
- channel of 850mm or shorter

Conventional Flexible Lumen Sterilization Cycle with Extra Loads • Ten(10) Polyethylene and Teflon

STANDARD MODE

- lumen instruments with: Internal diameter of 1 mm or larger and a length of 1,200 mm or shorter
- One(1) Long Polyethylene and Teflon lumen instruments with: Internal diameter of 1mm or larger and a length of 4,000mm or shorter

Versatile Flexible/Rigid Sterilization Cycle with Extra Loads

- Ten(10) Stainless steel lumen instruments with: Internal diameter of 1 mm or larger and a length of 500 mm or shorter
- Five(5) Stainless steel lumen instruments with: Internal diameter of 0.77 mm or larger and a length of 500 mm or shorter

One shelf of the chamber can be used in SURFACE mode.

Max. loading weight: 7.4kg (16.3lbs)

Up to three shelves of the chamber may be used in ENDO mode.

*Only Flexible Endoscopes with accessories may be loaded in ENDO mode

Up to two shelves of the chamber may be used in STANDARD mode.

Max. loading weight: 14.8kg (32.6lbs)

Up to two shelves of the chamber may be used in ADVANCED mode.

Max. loading weight: 14.8kg (32.6lbs)

No.of cycle/H2O2 bottle

No.of cycle/H2O2 bottle

No.of cycle/H2O2 bottle

No.of cycle/H2O2 bottle

^{*} All medical devices to be sterilized must be washed and dried completely prior to sterilization in HMTS system.



HMTS-142SE

HMTS-142 SE is the latest and most advanced Low Temperature Hydrogen Peroxide Sterilizer from Human Meditek. Designed to fulfill the evolving global customer demands, HMTS-142 SE offers unique features to make Central Sterile Supply Departments more productive, dependable and safe.

ıııı Efficient ıııı

With four cycle modes, HMTS-142 SE can sterilize a wide variety of medical instruments quickly and efficiently. HMTS-142 SE features an enhanced dual plasma sterilization technology powerful enough to achieve sterility assurance level (SAL) of 10–6 at half cycle while the large capacity 142L chamber allows you to process more loads per cycle.

um **Safe** um

HMTS-142 SE sterilizes with lower concentration of hydrogen peroxide (50%) which reduces the risk of damage to expensive medical devices. The chemistry safely breaks down into water and oxygen without leaving any toxic residues, ensuring optimal safety for users and the environment. In addition, HMTS-142 SE is equipped with Independent Monitoring system (IMS) which guarantees sterilization parameters are met every cycle.

nm Easy to use nm

HMTS-142 SE features aesthetically pleasing design with a large display screen that is easy to read. Enhanced user interface and graphic loading instruction reduce the potential for human error. Ergonomically placed automatic sliding door with foot pedal allows hand-free loading and the L.E.D indicator lets the users to monitor the sterilization cycle even from a distance.



Colorized L.E.D indicator



Intuitive and straightforward user interfaces





* Launch Scheduled at the end of 2018.

Accessories for HMTS-Series Low Temperature Hydrogen Peroxide Sterilizer

Consumable products and accessories are specially designed to ensure the highest compatibility and maximum validation efficacy for HMTS-Series Low Temperature Hydrogen Peroxide Sterilizers.

Sterilant Agent

HMSA-80 Sterilant, specially designed for HMTS-Series, contains 50% Hydrogen Peroxide. The lower concentration of hydrogen peroxide reduces the device and intrument damages while providing maximum sterilization efficacy. The self-contained bottles ensure user safety and can provide up to 18 cycles between changes.



Biological Indicators

HMBI-T is a self-contained biological indicator innoculated with Geobacilus stearothermophilus, a highly resistent bacterial spores, to safely and effectively monitor the sterilization processes.



Chemical Indicator (Strips / Tapes)

Designed specially for HMTS-Series Low Temperature Hydrogen Peroxide Sterilizers, The chemical inditors offer a distinct and easy to read color changes to provide confidence for proper sterilization.



Tyvek® Rolls

We offer a wide variety of sizes and styles of Tyvek® Rolls, which provides excellent tear and puncture resistance. The pouches and rolls can provide sterile storage for up to one year.



Other Accessories

We offer a variety of exclusive accessories, specially designed for HMTS-Series Low Temperature Sterilization System, to enhance convenience and efficiency.



*Human Meditek recommends the use of approved consumables and accessories to fully gurantee the sterilization processes.

HMTS has the Best Answer!

